

Response

Response


The Response is a versatile, accurate, hygienic and easy to use bench-top

revolutionary head designed to handle liquids, creams, pastes, and solids in

hoppers and other accessories, the

Response can be made to meet your

These machines are ideally suited to laboratory, R&D and small batch



Output (per minute)	< 40	< 80	Dependent
Depth			

Response

Compact and robust in a stainless steel case with clear, simple



Tube seal range	10–50 mm diameter
Output (per minute)	Up to 4
Electricity	
	400 × 345 × 700 mm
	A: fA™Y1BEA™yBØB1B'2BÄE Y™Ä B' 2a ØB1B'Ef B
	EBØ™Y%B‡ BÄ' {É' ' { Éa' AB9
	{BØØBÄÉ'' %Ø' %B‡ 1BÄÉ' { Éa' AB9

2 Y22Ø 2

Designed to close collapsible metal tubes, the Centrac FCC is operated by a handle, which requires

tubes in the correct position, a gauge plate regulates

3101000	
Output (per minute)	Up to 10
Tube diameters	
3110000	



Ø31BY™ 9P&Z Ä(rBÄE	
Output (per minute)	Up to 10
Tube diameters	203 mm

Ø Ä · Ä

